



Support Stand with Rod

OVERVIEW

The support stand with stainless steel rod accommodates Branson's converter and accessories. The converter clamp supports the ultrasonic stack (converter and horn) and is easily adjusted to properly position the horn in the sample.

The stand features a black enameled cast-iron base with zinc-plated rod. A full acoustic sound enclosure is also available for sound-sensitive environments.

STAND MEASUREMENTS

Model 8800

Stand Base: 9" l x 5.5" w x 24" h (229 x 140 mm)

Rod: 24" h x ½" dia. (609mm x 13mm)

Branson Part No. 101-063-1108
Stand with rod

Branson Part No. 101-063-1108
Converter clamp for use with Branson SFX150, SLP Sonifiers

Branson Part No. 101-063-1110
Converter clamp for use with Branson SFX250, SFX550, 250A, 250D, 450A, 450D Sonifiers



Americas

Branson Ultrasonics Corp.
41 Eagle Road
Danbury, CT 06813-1961
T: 203-796-0400
F: 203-796-9838
www.Emerson.com/Branson

Europe

Branson Ultraschall
Niederlassung der Emerson
Technologies GmbH & Co. OHG
Waldstrasse 53-55
63128 Dietzenbach, Germany
T: +49-6074-497-0
F: +49-6074-497-199
www.Emerson.com/Branson

Asia

Branson Ultrasonics (Shanghai) Co., Ltd.
758 Rong Le Dong Road
Song Jiang, Shanghai, PRC, 201613
T: 86-21-3781-0588
F: 86-21-5774-5100
www.Emerson.com/Branson